



LIMITED TIME PROMOTION THROUGH MAY 31, 2024.

Save now on our best-selling products...

Be sure to check all our promos. We've put together some special offers on the products you use the most.

EMS High Precision Ultra Fine Tweezers

These very popular, high-quality tweezers are ideal for laboratory work, including precise handling of TEM grids. They have an anti-glare satin finish. SA low carbon austenitic steel is non-magnetic and anti-acid. It is ideal for use where corrosion resistance and toughness are primary requirements.

EMS High Precision Ultra Fine Tweezers, Style 3X

Tip Configuration: Very sharp, self closing.

SPECIAL OFFER: BUY 2 UNITS, GET A THIRD FREE.

Description	Material	OAL	VWR Cat. No.	Unit
Style 3X	SA steel	120 mm	100494-728	Each

EMS High Precision Ultra Fine Tweezers, Style 4X

Tip Configuration: Extra fine, self closing

SPECIAL OFFER: BUY 2 UNITS, GET A THIRD FREE.

Description	Material	OAL	VWR Cat. No.	Unit
Style 4X	SA steel	110 mm	100494-738	Each

Dumont Tweezers, Style 5


- Hand made of the highest quality steels and alloys
- High precision Swiss manufacturing
- Perfectly finished to suit your application

Tip Configuration: Straight positive action with super thin tips.

Dumoxel® has excellent resistance to sulphuric environments, hydrochloric acid, as well as all other mineral and organic acids. Dumoxel® is 95% antimagnetic, resistant up to temperatures of around 400°C (DIN 50 914) and is suitable for autoclave sterilization at 270°C.

SPECIAL OFFER: BUY 2 UNITS, GET A THIRD FREE.

Description	Material	OAL	VWR Cat. No.	Unit
Style 5, High Precision	Dumoxel®, Polished	110 mm	76375-932	Each



Purchase the qualifying product(s) by visiting vwr.com, calling **1 800 932 5000**, or contacting your VWR Sales Representative between now and May 31, 2024. To redeem your free product, have your proof of purchase ready (VWR invoice or packing slip), visit vwr.com/promotions, and search for **Promo Code 6049**. Click "redeem" and complete the online form before June 30, 2024.