

# Tissue Homogenizer and Cell Lyser

Promo Code 6089

Receive a starter pack free with every CG-600 purchased.

Qualifying Equipment	
Product	VWR Cat. No.
HG-600 Geno/Grinder 2010 Tissue Homogenizer and Cell Lyser, 115 VAC, 60 Hz	76181-808

Choose ONE of the packages below and redeem on [vwr.com/promotions](http://vwr.com/promotions)

## Package 1 – High-Capacity Clamp and Ceramic Grinding Cylinder Starter Pack

Item #	Accessory Description
2014	<b>High-Capacity Clamp</b> – Thick, rugged foam block holds 30 standard 50 mL centrifuge tubes. Includes clamp insert (foam holder, bottom, lid and locking knob). Fits into the Geno/Grinder 2010 clamp (no assembly required)
2183	<b>Large Ceramic Grinding Cylinder</b> – Ceramic grinding cylinder with angle-cut ends. Package of 100. Chemically inert ceramic cylinders. The angle-cut ends provide better homogenization in applications such as QuEChERS pesticide residue extraction, reducing time and resulting in a thorough extraction of analytes into the solvent.

## Package 2 – Jar Starter Pack

Item #	Accessory Description
2259	<b>Holder for PET Jars</b> – Foam holder for two 25 oz. (750 mL) polyethylene terephthalate jars (#2258) or four 12 oz (355 mL) polyethylene terephthalate jars (#2248)
2258	<b>Jars</b> – 25 oz. (750 mL) polyethylene terephthalate (PET) jar. Package of 12. Holds ten #2156 11mm stainless steel grinding balls. Grinding load 60 to 80 g sample per jar
2248	<b>Jars</b> – 12 oz. (355 mL) polyethylene terephthalate (PET) jar. Package of 16. Holds five #2156 11mm stainless steel grinding balls. Grinding load 15 to 25 g sample per jar
2156	Grinding Balls, nominal diameter 7/16 in. (11 mm) - Made of 440C stainless steel. Package of 100

## Package 3 – Titer Plate Starter Pack

Item #	Accessory Description
2205-10	96 Square Well Reinforced Titer Plate, Pack of 10
2206-10	Cap Mat for 2205, Pack of 10
2012T	Nesting Tray
2100	Grinding Ball Dispenser – Simultaneously dispenses one 5/32 in. (4 mm) steel grinding ball (2150) into each well of a standard 96-well titer plate
2150	4mm Stainless Steel Grinding Balls, 5/32 in., bag of 500.

Purchase by calling 1.800.932.5000, visiting [vwr.com](http://vwr.com) or contacting your Avantor sales representative. Have your VWR purchase order, VWR invoice or VWR packing slip ready and fill out the online form at [vwr.com/promotions](http://vwr.com/promotions) and search for Promo Code 6089 before October 31, 2024. Offer valid on orders placed between January 1, 2024 and September 30, 2024. Redeem before October 31, 2024. While supplies last. Allow up to 8 weeks to receive your free product. Void where prohibited.