

Vacuum Pump Trade-up

VACUUBRAND® Vacuum Pumps

It's time for an upgrade! Take advantage of special pricing through VWR and SAVE on VACUUBRAND products.

VACUUBRAND oil free diaphragm pumps offer convenience, performance, reliability and value for your laboratory vacuum needs. Trade-up offer applies to select pumps (listed below).

To qualify for special Trade-up pricing, place an order with the VWR channel with any of the part numbers listed below, then e-mail promotions@brandtech.com, referencing the avantor Trade up promotion and indicate what vacuum pump you replaced.

VACUUBRAND Oil-Free Diaphragm Vacuum Pumps

- Chemistry design diaphragm pumps incorporate chemically resistant wetted materials for excellent corrosion resistance against aggressive gases and vapors.
- Extremely low maintenance. There's no oil to change or monitor! Typical diaphragm service life up to 15,000 hours of use—that's years in most applications.
- Well suited for most laboratory applications including, rotary evaporators, centrifugal concentrators, vacuum ovens, solid phase extraction and vacuum filtration.
- VACUUBRAND diaphragm pumps are whisper quiet. Most laboratory models operate at about the sound level of a quiet conversation.



BT TradeUP Cat. No.	VWR TradeUP Cat. No.	Applications	Description ("C" denotes Chemistry design)	TradeUP Price
20721003TU	MP97058-164	Filtration, SPE	ME 1, 100 mbar, 14 lpm	\$890.00
20721103TU	MP97058-166	Filtration, SPE	ME 1C, 100 mbar, 14 lpm	1,190.00
20724103TU	MP89412-354	Rotary Evaporation	MZ 1C, 12 mbar, 15 lpm	2,190.00
20732303TU	MP97007-180	Rotary Evaporation, Vacuum Ovens & Concentrators	MZ 2C NT, 7 mbar, 38 lpm	2,595.00
20696603TU	MP76529-726	Rotary Evaporation, Vacuum Ovens & Concentrators	MD 1C, 2 mbar, 25 lpm	3,495.00

*US plug included, all pumps work at 120V/60Hz.

Offer valid through 12/31/24

Promotional part numbers are available for a limited time, must be used at the time of purchase, and may be ordered via the "ORDER ENTRY" link on vwr.com. (Promotional part numbers cannot be found via Search on vwr.com.) For assistance with ordering, please contact Avantor Customer Service at 1.800.932.5000. All prices in US dollars.


avantor™
delivered by **VWR™**

VWR.COM

Prices, product, and/or services details are current when published and subject to change without notice. | Certain products or services may be limited by federal, state, provincial, or local regulations. | VWR, part of Avantor, makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC and/or Avantor, Inc. or affiliates. All prices are in US dollars unless otherwise noted. Offers valid in US, void where prohibited by law or company policy, while supplies last. | Trademarks are owned by Avantor, Inc. or its affiliates, unless otherwise noted. | Visit vwr.com to view our privacy policy, trademark owners, and additional disclaimers. © 2024 Avantor, Inc. All rights reserved.